

### Annex 3: analysis of breast milk samples in the laboratory

This appendix describes the procedures for the preparation and PCR testing of breast milk samples. All analyses were carried out at the laboratory of the National Institute of Public Health (INSP).

#### Sample preparation


Milk should be pre-processed to reduce its lipid and cellular component content.


Centrifuge the milk at 2500g (3500rpm) for 15 min

Recovering the supernatant

Dilute the sample (supernatant) to 1/10

#### Extraction





## Amplification

### Preparation of the reaction medium (Master Mix)

Briefly vortexing the reagents

Prepare the Master mix in the following proportions

Calculate the amount of PCR mixture as follows: Volume/ reaction\*N (samples + PC+NC+NTC+1)

Number of reactions	1	10
2X RADI FAST MasterMix	10µl	100µl
Mpox Primer and Probe Mixture	5µl	50µl
Volume du Master Mix	15µl	150µl
Samle, PC or NC	5	-

### Preparing the Reaction

Pipete 15µl of Master Mix into each reaction well of the PCR tubes

Add 5µl of sample (DNA extracts) or 5µl of control (CP, CN, or NTC)

Make sure all 3 controls are used for each handling

### Thermal Cycler Programming (CFX96)

This programming may vary depending on the type of instrument to be used

### Settings

Define the following

Settings	
Reaction Volume	20µl
Ramp rate	Default
Passive reference	None

### Selecting fluorophores

Target	Detector Name	Reporter
Monkeypox virus specific DNA Clade 2 (W. African)	MPXV	FAM
Clade1 (Congo Basin)	MPXV	Hex
IC	Internal Control	Cy5

### Temperature profile and acquisition

	Stage	Cycle repeats	Acquisition	Temperature [°C]	Time
Denaturation	Hold	1	-	95	20sec
Amplification	Cycling	45	-	95	2sec
			Yes	60	5sec

### Data analytics

For basic information regarding data analysis on a specific RT-PCR instrument, refer to the instrument's user manual.

### Interpretation of qualitative results

Sensing Channels			Interpretation of the results
FAM (Clade2)	HEX(Clade1)	Cy5 (IPC)	
≤40	>40 or -	≤40	Clade 2 Positive
≤40	>40 or -	-*	Clade 2 Positive
>40 or -	≤40	≤40	Clade 1 Positive
>40 or -	≤40	-*	Clade 1 Positive
>40 or -	>40 or -	≤40	Negative
>40 or -	>40 or -	-	**invalid. PCR inhibition or reagent problem. Resume the analysis from the primary sample or collect and analyse a new sample.
*: When the CPI Ct value is 40, the test is valid. In positive cases, IPC may not be detected due to competitive reaction with target genes. **: Results are invalid. A new test is necessary.			

### Quality control procedures

The quality of the sample is assessed during amplification. The amplification reagent contains a Sample control (SC), detected by the Cy5 fluorophore. In addition, the IC, always added at the beginning of the process (during extraction), allows you to control all the steps from extraction to amplification.

